



# YOUR ULTIMATE MILLING PARTNER


- 10 YEARS OF EXCELLENCE. A BOLD NEW BEGINNING. -



**2025 - 2026 GUIDEBOOK**

Service Price List

# HOIL DENTAL



HOIL DENTAL'S NEW  
HEADQUARTERS (HQ)

## CONTENTS

<b>04</b>	<b>Celebrating 10 Years of Excellence</b> Why Hoil? CEO Letter	<b>12</b>	<b>R&amp;D</b> 2024 Development 2025 New Services
<b>06</b>	<b>Introduction to Hoil Dental's New Journey</b>	<b>14</b>	<b>2025 New Services</b> Hoil Master Library Hoil Premium Bar Milling Service Hoil Smart Abutment SLM Partial Cobalt-Chrome Frame Hoil Multi Scanbody Hoil Scanbody
<b>08</b>	<b>2025 Price Freeze</b> Why Smart Labs Choose Hoil?	<b>26</b>	<b>Social Responsibilities</b> ESG Management Paperless Invoice, CIG Form
<b>10</b>	<b>Hoil Dental New Headquarters</b> Hoil Education Centre	<b>34</b>	<b>2025 Service Price List</b> Milling / Printing Services Available Milling Implant System
		<b>52</b>	<b>Consumables</b> Hoil Dental Shop

# CELEBRATING 10 YEARS OF EXCELLENCE



Since our establishment, Hoil Dental has been committed to innovation and excellence in the dental industry. As we celebrate our **10th anniversary in 2025**, we take pride in our position as the leading milling centre in the UK.


With **over 420,000 implant interface millings successfully completed** — all with zero connection failures — our expertise and precision are backed by a decade of consistent performance. In total, we have **completed over 3 million cases** across a wide range of specialties, including implants, zirconia, 3D printing, and more. This milestone reflects not only our technical capability but also the trust and support of our valued partners.

In 2025, we have invested heavily in new technology to further enhance our capabilities. This includes the launch of the **Hoil Premium Bar Milling Service, Hoil Smart Abutment, SLM Partial Cobalt-Chrome Frame, Hoil Multi Scanbody, Hoil Scanbody**, and the **opening of our new headquarters featuring an education centre**.

A wide range of lectures and training programs will also be introduced as part of this initiative, strengthening our commitment to quality, innovation, and professional growth. We aim to empower our partners with knowledge and technical excellence - reinforcing our shared vision for the future.

## WHY HOIL?

*"10 years of **dedication, innovation,**  
**and impact in Hoil**"*



## CEO LETTER



**Hoil Choi**  
Chief Executive Officer

### Hoil Dental Celebrates 10th Anniversary with New Office and Cutting-Edge Technology to Deliver Premium, Affordable Care

As we celebrate our 10th anniversary this year, I am filled with immense pride and gratitude for the journey we've taken together. What started as a small, passionate team has now grown into a leader in the dental industry, all thanks to your trust and support.

At the heart of our vision is a commitment to making exceptional dental care accessible at a reasonable price. We believe that everyone deserves top-tier service without compromise, and it is our mission to ensure that you receive the very best, both in terms of quality and affordability.

This milestone is not just about looking back at our achievements, but also about looking forward to the exciting future we are building. This year, we have relocated to a new, state-of-the-art office, a space that reflects our commitment to innovation and providing the best possible experience for our customers.

Thank you for being a part of our journey. We are excited about the future and look forward to continuing to serve you with excellence for many years to come.

Sincerely,

Over the past decade, we have made a substantial investment in cutting-edge technologies, driving the development of advanced solutions. This investment allows us to offer premium services that were once considered out of reach, ensuring that we deliver not only the highest quality dental care but also unmatched value to our clients.

Hoil Choi



# Where Milling Meets Implant-Level Precision



*The precision machines **trusted in implant manufacturing**—  
now **shaping the future of dental milling**!*

At Hoil Dental, precision defines our craft. **Leveraging Swiss-type lathe technology** - renowned for ultra-fine tolerances and impeccable finishing - we elevate the standard of dental milling, particularly in our advanced applications. This means our customers get **implant-grade quality** in every part, every pass, every time.

We deliver precision you can trust, even with the most complex dental prosthetics.  
We believe that our growth depends on our customers. With sustainable innovation and honest value, we promise to deliver the best. Together with our partners, let's build a future together, shaped by precision and trust.



## STEPPING FORWARD AS AN ORIGINAL IMPLANT MILLING CENTRE

### Introduction to Hoil Dental's New Journey

#### *A Decade, a Million Milestones, and a **New Leap Forward***

Over the past 10 years, Hoil Dental has specialised in the precision processing of a wide range of original implant systems and compatible components, building an extensive repository of technical data, clinical know-how, and customer insights.


Grounded in this wealth of expertise, we are proud to have been officially recognised as an authorised milling centre through our collaboration with Arum Implant; a solution that reflects our dedication to precision, reliability, and innovation.






# 2025 PRICE FREEZE EVEN BETTER VALUE FOR LABS

## 2025 PRICE INCREASE COMPARISON



In a time when inflation is forcing many to increase prices, we've taken a different path. Through relentless efforts in **R&D and production automation**, we've managed to lower our costs — so we can continue offering you the same prices, with even better quality.



## INVEST IN WORK EFFICIENCY

### & FOCUS ON YOUR LAB GROWTH

*Buying a milling machine is easy,  
but **running one is a whole different story.***

Category	In-House Milling	Hoil Milling Centre
Metal Equipment Cost	£75,000 + VAT (Average)	£0
Maintenance & Tools	£10,000 (Average per year)	£0
Skilled Technicians	£45,000 / per person (at least one)	Not Needed
Software Update & Annual fee	£1,800 (Average per year)	£0
Library Update	£500 - £700 / per implant	£0
Focus on Design	Limited (divided focus)	100% Focused
Technology & Skill	Restricted system Quality stagnation without reinvestment	900+ System Improved quality via reinvestment
Certification	X	ISO13485, UKCA, CE UK Exclusive angulated licensee
Gold Coated Screw (i-ON)	X	Free
Material Cost		



Custom Abutment

£50 - £60 + α

\* Avg. Premill + Screw price is £60.

£47.95



All-on-4 Bar

£246 ~ / per bar + α

\* This reflects the average cost when milling 20 bars per month.

£230

## WHY SMART LABS CHOOSE HOIL?




Material cost  
Stock management  
Hiring & training skilled technicians  
Maintenance, tools, extra costs  
Regular software updates  
Machine downtime

All Solved! by

**HOIL**  
DENTAL





#### Engineered for Innovation and Quality

Our newly upgraded production facility is equipped with state-of-the-art dental milling machines and industrial machines, with dedicated spaces for each team to monitor and oversee every step of the production process.

# HOIL DENTAL NEW HEADQUARTERS

At Hoil, we have always looked forward to connecting with our customers—not simply to deliver products, but to build true partnerships.

*“In what ways should we connect, and what should we share?”*



In response to this question, we have established **our new office as the Hoil Education Centre**—a dedicated space for education, collaboration, and innovation. This centre serves as the foundation for delivering meaningful learning experiences that go beyond traditional product support.

As part of this initiative, **we offer two distinct educational programmes:**

The first is an advanced masterclass, featuring globally recognised speakers who deliver high-level insights into the latest trends and techniques.

The second is an introductory digital dentistry course, specifically designed for those seeking to build a strong foundation in digital workflows and technologies.

Through these initiatives, we aim to foster growth, expertise, and innovation across the dental community.

## *A Partnership that Grows Together - Education Centre*

We're planning to host various lectures and hands-on courses regularly. The space is also available for free rental when needed—so feel free to reach out to us anytime if you'd like to use it!

[info@hoildental.com](mailto:info@hoildental.com)



*Follow Hoil Dental on social media for the latest updates!*





The world widest implant systems available  
**900+** by Hoil R&D team

Driven by our R&D team’s commitment to innovation, we have developed significantly more milling implant systems than last year, now offering **over 900 systems**.

We continuously expand our available milling implant systems to quickly and accurately respond to the diverse needs of our customers. With the launch of every new implant system in the market, we have developed optimal solutions to address customer demands. As a result of these efforts, our library now encompasses **more than 900 systems**, with **24 new systems added just last year**.

Our R&D team works closely with the CS team, carefully listening to customer inquiries regarding new implant systems. Based on this input, we consistently develop and expand our offerings. Building on the foundation of over 900 systems, our library is steadily moving toward the 1,000 mark. This number is more than just a figure—it represents the trust of our customers and our dedication. It is our pride.

You are welcome to request the **development of any new implant system for milling!**

Our R&D team is dedicated to meeting each and every client need, ensuring that no request goes unaddressed. Even if you need just a single case, **please feel free to reach out** — we are here to provide you with proactive solutions.

support@hoildental.com

2024 Development  
Compatible 19 Implant Systems, 24 Sizes

Mar	1	BTI multi 3.5
Apr	4	MIS® V3 Standard Sweden&Martina® PAD.r Multi-Unit STRAUMANN® BLX® 5.0/5.5/6.5(WB) H ZIMMER® Universal Abutment Multi-Unit
May	3	ADIN® UniFit™ CAMLOG® iSy® MIS® V3 Narrow
Jun	3	Biohorizons® Tissue 3.5 Biohorizons® Tissue 4.5 Biohorizons® Tissue 5.7
Aug	1	IHDEDENTAL Hexacone®
Sep	3	BIODENTA® Bone Level B1 3.5 BIODENTA® Bone Level B2 4.1/4.8/6.0 Osteo Ti Internal Hex
Oct	4	BEGO® Internal Hex 3.0 BIOTEC® Multi-Unit B1 3.5 BIOTEC® Internal Hex Tree-Oss® Internal Hex
Nov	3	Klockner® SK2 MEGAGEN® Blue Diamond NC MEGAGEN® Blue Diamond RC
Dec	2	ZIACOM® NP 3.2/3.3 MIS® Connect

2025 New Services



Hoil Master Library ..... Page. 14-15  
- An Industry Largest



Premium Bar Milling Service ..... Page. 16-17  
- Swiss-type Industrial Milling System



Hoil Smart Abutment ..... Page. 18-19  
- Direct Abutment



SLM Partial Cobalt-Chrome Frame ..... Page. 20-21  
- 3D Metal Printing



Hoil Multi Scanbody ..... Page. 22-23  
- Intraoral Bar Scanbody



Hoil Scanbody ..... Page. 24-25  
- Intraoral Single Scanbody



Launched Feb. 2025

# HOIL MASTER LIBRARY

- ✓ COMPATIBLE WITH ANY INTRAORAL/MODEL SCANBODY YOU USE
- ✓ ORDER AT HOIL MILLING CENTRE USING MASTER LIBRARY
- ✓ FULL USE OF HOIL SHOP COMPONENTS

## Universal Compatibility

Since the launch, our customers have reported increased satisfaction with the ability to integrate ARUM digital analogs with a wide range of intraoral and model scanbodies. This flexibility allows dental professionals to choose the tools that **best suit their practice without compatibility concerns.**



### Happy Customer

I've been waiting for this when I saw it on the website. Thank you HOIL!



### Happy Customer

Cool idea! I love HOIL! You're the best!!



### Happy Customer

We are happy to see this new library that frees us from intraoral scan data restrictions with multiple scanbody brands!



INCORPORATE DIGITAL WORKFLOW

TIME SAVING AND REDUCING ERRORS

A significant milestone in our mission to enhance the dental laboratory experience



**2 CLINICAL  
SCREWS  
INCLUDED**  
**STRAIGHT**  
Same as the original  
screw design  
**ANGLED**

# COMPATIBLE WITH MAJOR IMPLANT INTRAORAL / MODEL SCANBODY

# HOIL PREMIUM BAR MILLING SERVICE

*Invested in **Technology**, Committed to **Perfection**.*

Every clinical case is unique, and so are our customers' designs. Traditional milling systems often struggled with complex STL geometries due to undercuts and machining limits. To break through these limitations, we made a substantial investment in the advanced **Swiss-type industrial milling system**—engineered to handle even the most intricate designs with precision.

We no longer face limitations when it comes to **special design production**. Whatever our customers envision, we can now deliver **accurately, consistently and without compromise**.



**TRANSFORMING DENTAL COMPONENTS  
INTO PRECISION ART**

*We Don't Just Mill Parts  
— **We Engineer Innovation***



## ✓ Ultra-Precision & High-Speed 5-Axis Machining

*Swiss-made high-precision machine.*

*Heidenhain TNC 640 multi-touch controller for intuitive operation.*

*High-Speed Machining (HSM) technology for superior surface quality and fit*

## ✓ Exceptional Accuracy & Rigidity


## ✓ Optimised for complex geometries, custom dental bars and frameworks

## ✓ A wider range of angulated milling solutions

## ✓ Precision milling + measurement processes ensure flawless results

## ✓ Ensure ultimate patient comfort and functional precision


***Ultimate Solutions**  
**Tailored to Your Needs - Page. 37***



*Premium Telescopic Bar*



*Premium Basic Bar*



*Premium Hybrid Bar*



*Premium Abutment Hybrid Bar*



# WHAT IS DIFFERENT ABOUT HOIL SMART ABUTMENT?

(Direct Abutment)

Crafted from premium-grade titanium, the Hoil Smart Abutment sets a new standard in precision and prosthetic accuracy.

Unlike conventional custom abutments which are produced using premilled blanks, our abutments are **milled directly from solid titanium bars using high-end Swiss-type lathes**.

This advanced direct milling process delivers exceptional fit, flawless surface finish, and eliminates the need for additional rescanning.

## 1 100% in-house production ensures full quality control

Bar-based precision cutting allows for high reproducibility across batches, maintaining surface roughness and optimal geometry. Our customers can rely on predictable, uniform outcomes with no need for adjustments.

## 2 Swiss-type ultra-precision Lathe machining

Trusted by leading implant manufacturers worldwide, this machine delivers aerospace and medical-grade precision machining.



# Head-to-Toe Technology

without *Pre-milled Blanks*



## 3 No rescanning Required



Our precision bar-based manufacturing ensures exceptional fit accuracy with final crowns from the first scan. It allows direct milling without requiring a rescanning process, which simplifies the workflow and saves time. Therefore, the best adaptation quality can be pursued even if the crown and abutment are performed simultaneously in CAD.

## 4 Superior Margin Line


With our Hoil Smart Abutment, fine anatomical margins are reproduced more smoothly and accurately, improving soft tissue response and long-term prosthetic success.

## 5 Better Precision and Fit

Hoil Smart Abutments are milled entirely from a solid titanium rod in a single CAM process. This allows for more accurate control of geometry, resulting in a better fit and reduced risk of misfit between the abutment and implant.







## DEDICATED DENTAL METAL 3D PRINTING

High-end machine with advanced algorithm and components brings the fastest printing speed


## 3D METAL PRINTING

### SLM PARTIAL COBALT-CHROME FRAME

Once again, we, Hoil Dental, have made a **bold investment in 3D metal printing** to better serve the partial cobalt-chrome frame customers in the UK. We have committed a substantial investment to strengthen our technological leadership in the metal printing sector.

- ✓ *Outstanding **Speed and Efficiency***
- ✓ *Exceptional **Cost Effectiveness***
- ✓ *Proven **Flexibility and Adaptability***

## WORKFLOW



### PARTIAL CO-CR FRAME PACKAGE SERVICE

- |           |                                      |                  |                  |
|-----------|--------------------------------------|------------------|------------------|
| Option 1. | Printing + Polishing                 | +2 Days Dispatch |                  |
| Option 2. | Scan / Design + Printing             | +2 Days Dispatch |                  |
| Option 3. | Scan / Design + Printing + Polishing | £10 off          | +3 Days Dispatch |



# HOIL MULTI SCANBODY

**NEW!**



✓ OPTIMISED ACCURACY OF INTRAORAL SCANS FOR FULL-ARCH SCANNING

✓ COMPATIBLE WITH ALL INTRAORAL SCANNERS

▲ **Direct CAD Workflow Support**

Direct design in Exocad and 3Shape without separate scan merging

◆ **Seamless scan continuity**

ensured by minimising gaps between scanbodies

◆ **Enhanced positional alignment**

with a multi-faceted polygonal design

■ **Smooth, rounded edges**

for optimal scanning and a comfortable intraoral experience

● **Biocompatible materials**


gold-coated titanium for superior safety and durability

▶ **Easy to Clean**

Autoclavable for safe reuse, without compromising scanning accuracy even after multiple sterilisations

◆ **Efficient System Identification**

- High-precision laser markings on the scanbody enable rapid and accurate identification of the corresponding implant system
- Designed for seamless integration with multiple implant systems, ensuring broad compatibility and workflow efficiency



COMFORTABLE EXPERIENCE

TIME SAVING AND REDUCING ERRORS

INCORPORATE DIGITAL WORKFLOW



**I N T R A O R A L**  
**POSITIONING SCANBODY FURTHER**  
**INWARD OPTIMISES SCAN ACCURACY**

PRECISION

SIMPLICITY

DIGITAL POWER



# I N T R A O R A L

## WHERE PRECISION MEETS DIGITAL SIMPLICITY AND ACCURACY

**NEW!**

## HOIL SCANBODY

### *The Optimal Choice for **Digital Implant Workflow***

In the rapidly evolving world of digital dentistry, clinicians increasingly demand **faster workflows, greater accuracy, and seamless integration with various implant systems**. Recognising these challenges from direct feedback by our valued customers, we set out to engineer a solution that met the real needs of the modern dental practice.

*"We didn't just want to make another scanbody — we wanted to create a tool that truly simplifies digital implant workflows without sacrificing accuracy."*

- R&D Team, Hoil Dental -



#### **Reflection-Free Scanning**

- Special matte surface finish eliminates light reflection for clearer scan data
- Enhances scan accuracy and eliminates the need to apply spray to the scanbody



#### **Convenient System Identification**

- Laser markings on the scan body enable easy and fast system identification



#### **Optimised Anatomical Design**

- Minimises diameter at the gingival and alveolar bone regions for versatile case application
- Provides excellent adaptability based on implant placement depth



#### **Wide Intraoral Scanner Compatibility**

- Compatible with all major intraoral scanners
- Optimised for model-less digital workflows



#### **Material and Durability**

- Manufactured from high-purity titanium for outstanding biocompatibility
- Maintains structural integrity for long-lasting use, even after repeated sterilisation



#### **User Convenience and Stable Fixation**

- Screw-retained fixation provides a secure and stable connection to the implant fixture
- Offers easy removal after scanning, enhancing clinical workflow efficiency



#### **High Precision and Accuracy**

- Precision-engineered design ensures accurate scanning of implant positions and angles
- Maximises consistency and reliability throughout the CAD/CAM production process



#### **Optimised for Digital Workflow**

- Enables immediate CAD design after fast and accurate scanning
- Reduces production time and improves quality compared to traditional analog methods



# Lead the way in Sustainable Business

Commitment to innovation while maintaining a strong focus on sustainability

Our dedication to sustainable practices is part of a broader corporate ethos that emphasises delivering value through cost-effective solutions. By reducing operational costs associated with paper usage, we can pass these savings on to our customers, making quality dental solutions more accessible.

Powered by 100% Zero Carbon Footprint — Certified and Verified

All electricity used at Hoil Dental is fully backed by renewable energy and nuclear power sources, certified through Renewable Energy Guarantees of Origin (REGOs) and Nuclear Declarations. This ensures that our energy supply is entirely supported

by zero carbon electricity, enabling us to operate with no associated carbon emissions. Through this commitment, Hoil Dental proudly advances its ESG vision, contributing to a cleaner, more sustainable future.

CERTIFIED  
100%  
ZERO CARBON  
FOOTPRINT

28.8  
Trees saved

600,000L  
Water saved


1.2 TONS  
CO<sub>2</sub> reduced

# WE VALUE THE ENVIRONMENT

ESG | Environmental Social and Governance

PAPERLESS INVOICE / SOLAR PANEL / METAL SLUG RECYCLING

# SWITCHED TO DIGITAL DOCUMENTS FOR A SMARTER, GREENER FUTURE



## Paperless invoice, CIG form


Launched Jan. 2025

We have implemented **electronic invoices** and a **digital Case Instruction Guide (CIG) form**, effectively minimising paper use and promoting eco-friendly solutions. This move not only reduces the environmental impact of their operations but also enhances **efficiency in communication and workflow** management. By transitioning to digital formats, we aim to streamline processes, allowing for quicker responses and better organisation. This shift reflects the company's commitment to innovation while maintaining a strong focus on sustainability.

✓ Invoice & CIG NEW

Download on Hoil Order System

**CIG** | Case Instruction Guide



Tooth	Compatible Implant System & Size	Screw Channel	Item Code	Driver Size	Torque	Special Technical Advice for User
16	Sweden&Martina® Prama 4.25 NH	Straight	S668-G	HEX 1.3	25	
12	DENTSPLY® Ankylos® PILAR REGULAR 3.5 NH	Angulated	DS013-G	Ball Screw Driver Torx	15-20	Please adhere to the driver size and torque specifications. This angled screw is developed by Hoil Dental Milling Centre.
22	DENTSPLY® Ankylos® PILAR REGULAR 3.5 NH	Straight	S136-G	HEX 1.0	15	

### ✓ Clear Instructions for Each Tooth Unit

Especially when different implant systems are used depending on the position, the CIG specifies exactly which screw is required for each tooth unit. This helps avoid any mix-ups and ensures clarity for both the lab and the dentist.

### ✓ Driver Size & Torque Guidelines

The CIG includes important notes on driver size and torque. Please make sure to review these instructions carefully when using the screw to ensure proper handling and avoid complications.

### ✓ Lifetime Warranty Documentation

The CIG can also serve as a formal case record for dentists, providing a reliable document for lifetime warranty claims.

# WORK PROCESS

STL File Entry Service Lead Time

- ✓ *Fast Turnaround*
- ✓ *Quality Work*
- ✓ *Technical Support*



## Same day Dispatch

Order by 10AM

SLM 3D Metal Printing



Order by 11AM

Metal Milling



Zirconia Prime Service (up to 3-unit bridge)



## Next day Dispatch

Order by 10AM

SLM Partial Cobalt-Chrome Frame



i-ON Coating Service



Order by 12AM


Zirconia Standard Service





# DESIGN & MILLING SERVICE

Next Day Dispatch Timeline



## DAY 1


### 1. Sending us your model

Send us your model with an order form **before 5 PM**. If it is over 6-unit, please send model with denture setup or wax-up.



### 2. Model scan

Scan the working model by using ARUM 3D scanner and ARUM scanbody. This stage is very important in order to get the best results.




### 3. Designing

Design the dental restoration digitally with CAD software by experienced dental technicians.



### 4. Sending you design files


After designing, we will send design files of your job for confirmation the next day **before 10 AM**.



## DAY 2

### 5. Design confirmation

Confirm our design approval email before milling stage. Please accept it or notify us anything further **by 11 AM**. If not, dispatch date will be delayed.



### 6. Milling

Once it is confirmed, we will mill the restoration using our state of the milling machines and Metal, Zirconia, PMMA, PEEK, Composite, etc., achieving high fidelity.



### 7. QC Check

Our experienced senior dental technicians check the final quality control.



### 8. Packing & Shipping

Pack your jobs carefully with an order slip, then ship to you with a trackable express courier service.






# AI POWERED AUTO ENTRY SYSTEM

We developed Auto Entry to eliminate the repetitive, time-consuming process of manual data entry. By applying AI technology, **our goal was simple: to help you save time, reduce errors, and work more efficiently.**

With Auto Entry, our customers can complete tasks twice as fast, with greater accuracy and consistency. This not only improves productivity, but also helps reduce unnecessary issues caused by input mistakes.

## 4 TIMES FASTER, 100% SMARTER

*How much time can you  
save with **Auto Entry**?*



*Simply Drag & Drop Files*

Custom Abutment:  
Product: Titanium Custom Abutment  
Order Reference:  
'Click' or 'Drag and Drop' files here to upload  
Auto-Entry Required Files:  
exocad - stl, constructioninfo  
3shape - stl, xml  
☐ I hereby certify that the above statements are true and correct to the best of my knowledge.  
Confirm Order: ✓

1

### Save Time, Work Smarter

Manual data entry used to take up to two minutes per case. With Auto Entry, it takes just 30 seconds — cutting your workload more than half. Less time on repetitive tasks means more time for what truly matters.

2

### Increase Accuracy

Fewer typos. Fewer mistakes. Fewer corrections. Auto Entry dramatically improves data accuracy, helping prevent errors that lead to wasted time or rework.

3

### Smarter System, Constantly Evolving

Auto Entry is not an off-the-shelf tool — it's our own AI-powered system, developed in-house and constantly improved based on real customer feedback. That means you can always expect smarter, faster, and more seamless performance with every update.

## i-ON COATING


### Benefits for Hoil Milling Customers — Enjoy Hoil Dental's Exclusive Free Coating Service!

Since launching the service last year, we have seen a steady increase in demand for i-ON Coating. In line with our commitment to serving our customers, we continue to offer **FREE Titanium Nitride Coating exclusively for bar orders**. This allows patients to experience numerous aesthetic benefits. We encourage our valued customers to take full advantage of this exceptional service.


Additionally, we hope to widely inform dentists about the **FREE i-ON Screws** provided with all milling orders and that they can fully enjoy the benefits of coated screws.

- ✓ *Significantly Prevents Screw Loosening*
- ✓ *Enhanced Antibacterial Properties*
- ✓ *Improved Durability*

COATING SERVICE USER RATE



*Free i-ON Screw provided  
with Hoil milling*



*Free coating service for  
Hoil bar*

**Nanohardness at least 5 times,  
Elasticity at least 3 times more than uncoated!**

# 2025 SERVICE PRICE LIST



## Metal Materials

### Hoil Premium Medical Nitride Coating (i-ON Coating)

Titanium Nitride Coating technology serves to prevent framework contamination, significantly boosting **surface strength by 5 times** and extending **fatigue life by 3 times** compared to uncoated surfaces, as well as enhancing the products with an excellent quality finish.



Case	Price
Single Abutment (Ti, Co-Cr)	
Implant Bar (Anatomic bridge, Hybrid bar, Snap-in bar, Basic bar)	
Split Bar (Blender Bar)	

### Milling Service

### STL File Entry



#### Titanium Custom Abutment Inc. Screws



*Check Offered System*

offered system	+	Hex	+	Hoil/ Arum Scanbody	+	Straight
✓		✓		✓		✓

Non-engaging *or*  
Non Hoil/ Arum Scanbody *or*  
Non 'offered' implant systems

Angled Custom Abutment



#### Co-Cr Screw Retained Abutment Inc. Screws



#### All on Hoil - Multi Convertible Abutment (Each) Inc. 2 Screws + Thread Milling + Scan + Design

Opt. 1 - Free angle adjustment up to 50 degrees  
Opt. 2 - Adjustable gingival height & compression  
Opt. 3 - Adjustable multi-connection position

Shipping fee charges per parcel.





Hoil's Special Service - Bar Milling

- We charge the same price **regardless of the number of implants.**
- All bars come with a **free coating** upon request. (excl. Split Bar)
- **Re-coating services** are also provided **free of charge** if needed.

Milling Service

STL File Entry

 <div>2025 HIT</div>	<b>Split Bar (Blender Bar)</b> Inc. Screws	Titanium	2 - 8 implants
		Co-Cr	2 - 8 implants
	<b>Basic Implant Bar for denture</b> Inc. Screws	Titanium	2 - 8 implants
		Co-Cr	2 - 8 implants
	<b>Anatomic Bridge</b> Anatomically designed framework for composite, ceramic build-up Inc. Screws	Titanium	Pontic Position
			Implant Position
		Co-Cr	Pontic Position
			Implant Position
	<b>Hybrid Abutment Bar</b> Inc. Screws	Titanium	2 - 8 implants
		Co-Cr	2 - 8 implants

Shipping fee charges per parcel.






Hoil Premium Service

NEW LAUNCHING IN 2025



Milling Service

STL File Entry

	<b>Hoil Smart Abutment (Direct Abutment)</b> Inc. Screws	
	<b>Premium Basic Bar</b> Inc. Screws	
	<b>Premium Hybrid Bar</b> Inc. Screws	Full service only Scan + Design + Milling + Polishing
	<b>Premium Abutment Hybrid Bar</b> Inc. Screws	
	<b>Premium Telescopic Bar</b> Inc. Screws	Full service only Scan + Design + Milling + Polishing

Shipping fee charges per parcel.

Design Service Entry

We are offering design service from single unit to full arch. If it is more than 6 units, please send us the model with denture setup or wax-up. After designing, we will send you confirmation email.

Implant Case





Case	Scan + Design	Design
Implant position (per unit)		
Pontic position (per unit)		

Non-Implant Case

Case	Scan + Design	Design
Per unit		

Milling Service

STL File Entry

	Coping	Titanium	Single
			Bridge
	One-piece Locator Abutment Inc. Screws & LOC Retention Kit	Titanium	
		Co-Cr	
	One-piece Locator Bar Inc. Screws & LOC Retention Kit	Titanium	2 - 8 implants
		Co-Cr	2 - 8 implants
	Aluminium Bar for implant verification Jig: Up to 4 implant positions Inc. Screws		
			Additional implant position Inc. Screws

Shipping fee charges per parcel.

SURVEYOR  
HAND  
MILLING



Surveyor Hand Milling



Metal milling process only from model

Metal Bar + Zirconia Bar full package milling process from model

\* Please send us the model, this service is available when we start from scan and design.

Shipping fee charges per parcel.



# SLM 3D Metal Printing

## Partial Cobalt-Chrome Frame Service Entry

Printing	Scan + Design	Polishing	Price
✓			
✓	✓		
✓		✓	
✓	✓	✓	


\* Free coating service available upon request.

### Printing Service

### STL File Entry



Partial Cobalt-Chrome Frame



Co-Cr Coping  
Cement Retained

Single

Bridge  
up to 16 units



Co-Cr Full Crown  
Cement Retained

Single

Bridge up to  
16 units

Shipping fee charges per parcel.

## Design Service Entry

### Implant Case

Case	Scan + Design	Design
Implant position (per unit)		
Pontic position (per unit)		

### Non-Implant Case

Case	Scan + Design	Design
Per unit		

### Printing Service

### STL File Entry



Co-Cr Verification Jig  
Inc. Screws  
Only Try In

Up to 4 implants

5 - 8 implants



Co-Cr Secondary Structure



Co-Cr Dental Post

Shipping fee charges per parcel.

.....  
*Industry's First*

## One-day Trio Milling

*We value your time with the best quality!*



**Model (3D printed) + Custom Abutment + Zirconia milling**

“Experience faster turnaround times  
without compromising quality!”

5 times faster than ordinary sintering & Perfect sintering effect

Orders placed **by 11 AM,**

**Enjoy**  
**Same Day Dispatch!**



WhatsApp

## Soft Materials

### Milling Service

### STL File Entry

 Single or Up to 6 Unit Bridge Zirconia	Zr-I
	Supreme Symphony
	Hyperion
 7 or More Unit Bridge Zirconia	Supreme Symphony
	Hyperion
 Zirconia Hybrid Abutment Bar excl. Ti-base	3 implants
	4 implants
	5 implants
	6 implants
	7 implants
	8 implants
	9 implants

Shipping fee charges per parcel.



Design Service Entry

We are offering design service from single unit to full arch. If it is more than 6 units, please send us the model with denture setup or wax-up. After designing, we will send you confirmation email.

Implant Case






Case	Scan + Design	Design
Implant position (per unit)		
Pontic position (per unit)		

Non-Implant Case

Case	Scan + Design	Design
Per unit		

Milling Service

STL File Entry






	PEEK Hybrid Bar + Anatomically designed framework Inc. Ti-base	Pontic Position
		Implant Position
	PMMA	Multi-Layer PMMA Temporary Crown
		Multi-Layer PMMA Temporary Crown Bridge
	Composite Crown	Composite Crown
		Composite Crown Bridge
	Pink Gum Painting of the Zirconia (per unit)	
	PEEK Partial Frame Inc. Design	Upper
		Lower

Shipping fee charges per parcel.


3D Printing

Printing Service

STL File Entry

	Designing for Model Print	
	3D Printed Partial Arch	
	3D Printed Full Arch	
	3D Printed Gum per Unit	
	3D Printed Single Die per Unit	

Shipping fee charges per parcel.



# Additional Service

Printing Service

STL File Entry

Additional Service

STL File Entry

	Surgical Guide Partial Arch	Printing only	
		Inc. Design (without Metal Sleeve) 1 - 4 implants	
	Surgical Guide Full Arch	Printing only	
		Up to 4 implants Inc. Design (without Metal Sleeve) 5 - 9 implants	
	Surgical Guide Sleeves Prices may vary		

Shipping fee charges per parcel.


	Angulated Screw Channel for Metal Work	
	Angulated Screw Channel for Zirconia Work	
	Thread Milling for Bar	
	One-piece Locator (Customs Locator) Inc. LOC Retention Kit	
	Locator	LOC Bar Attachment
		LOC Retention Kit up to 20°

Shipping fee charges per parcel.



Additional Service

STL File Entry



Polishing



Laser Welding for Co-Cr and Titanium  
(Each welding)



i-ON Coating  
(Premium Medical Nitride Coating)



i-ON Screw

Shipping fee charges per parcel.

# HOIL MASTER LIBRARY

## THE WIDEST COMPATIBLE IMPLANT SYSTEM

*From A to Z, we've got every compatible implants system you need right at your fingertips!*



ADIN® CLOSEFIT™




ZIMMER®  
Universal Abutment Multi-Unit



*Become a Hoil Milling customer and  
unlock the Hoil Master Library now!*



Scan to Register!



# Available Milling Implant System

## A series compatible with

ADIN® CLOSEFIT™ / TOUAREG™ S&OS&SWELL / Multi-Unit® / UniFit™

ALPHA-BIO TEC® Internal Hex / CONCIAL / Multi-Unit®

ANTHOGYR® Axiom® BL / Axiom® TL / AnthoFit® OI / Multi-Unit®

ASTRA TECH™ OsseoSpeed™ TX / EV™ / EV™ PROFILE / Multi-Unit® / Angled / EV™ Multi-Base

AVINENT® OCEAN

## B series compatible with

BIO3 IMPLANT® Conical Connection / Multi-Unit®

BEGO® Internal / Multi-Unit®

BICON®

BIO COMP DENTAL INTERNAL HEX / Multi-Unit®

BIODENTA® Bone Level

BIOHORIZONS® External® / Internal® / Multi-Unit® / Tissue level

BIOMET 3i® External® / Certain® / Multi-Unit®

BIONER IMPLANT

BIOTEC® Internal Hex / Multi-Unit

BIOTECH® BIS 3.6 EXTERNAL / BIS Octagone Interne / Kontact / Multi-Unit®

BrainBase Corporation® MYTIS® Arrow Implant

BREDENT Medical SKY® / Copa SKY® / SKY® Multi-Unit®

BTI® INTERNA® / EXTERNA® / Multi-Unit®

BTK® Multi-Unit®

B&B DENTAL Conexa® / Multi-Unit®

## C series compatible with

CAMLOG® CAMLOG® / Multi / iSy®

CONOLOG® CONOLOG® / Multi

CHAMPION Multi

CLC conic™

CONEXAO® HI

CORE VENT Internal Hex

CORTEX® CORTEX® / Multi-Unit®

COWELLMEDI® INNO Internal

C-TECH® Bone level / Esthetic Line / Multi-Unit

## D series compatible with

DEEP®

DENTIS® s-clean / i-clean / Save / Multi-Unit®

DENTAURUM® TIO

DENTIUM® SuperLine / SimpleLine / NR LINE / Multi-Unit®

DENTSPLY® Ankylos® / Ankylos® Blance base / Xive® / Xive® Multi-Unit®

DIO® SM / UF / UF Multi-Unit

DYNA® PUSHIN Octalock®

## E series compatible with

EBI®

EDISON MEDICAL® Clicq Multi-Unit®

EUROTEKNIKA™ Aesthetica / Naturactis/Naturall+/Natea+ / Multi-Unit

## G series compatible with

GC Aadva™

Global D® Tekka® In-kone / In-Kone® Multi-Unit / EVL® Multi-Unit

GP Implant Multi-Unit

## H series compatible with

HIOSSSEN® ET / SS / Multi-Unit

HumanTech RATIO®

## I series compatible with

IHDEDENTAL Hexacone®

IPD Multi-Unit

IMPLANT DIRECT® Legacy™ / InterActive™ / Swish Active™ / Swish Plus™ / Tri-Lobe™ / Multi-Unit®

IMPLANT SWISS® Bone Level / Multi-Unit 4.8

INTRA LOCK GOLD & BLUE

iRES® Internal Hex 2.5mm / Internal Octa Tissue / Multi-Unit

## J series compatible with

J DENTAL CARE® JD Evolution® / JD Evolution® Plus+ / JD Icon® / JD Icon® Plus+ / JD OCTA® / Multi-Unit

We are updating the available implant list every month.

To get the most resent version of implant list, please visit our website, [www.hoidental.com](http://www.hoidental.com).

## K series compatible with

KENTEC®

KEYSTONE® Prima Connex®

Klockner® SK2

KYOCERA® FINESIA® Bone Level / Poiex

## L series compatible with

LASAK Bioniq® / Multi-Unit

LEADER Italia® Mores Connection

## M series compatible with

MEDENTiKA® Multi-Unit

Medentis Medical® ICX® Active Master / ICX® / ICX® Multi-Unit

MEGAGEN® Anyridge / AnyONE / Blue Diamond / Exfell / EZ Plus / MINI / MegaFix / RESCUE External / RESCUE Internal / Multi-Unit®

MIS® C1 / Connect / V3 / Internal Hexagon / Multi-Unit / Multi-Unit (New)

MSDI Alef-Spiral Internal Hex

## N series compatible with

NEOBIOTECH® IS SYSTEM / IT SYSTEM / IS ACTIVE SCRP

NEODENT® CM / GM / External / CM Conical / GM Conical

NEOSS Pro Active® / Access®

NOBEL BIO CARE® Active™ / Branemark® / Replace Select™ / Multi-Unit™ / Multi-Unit™ AESTHETIC/ N1™

NORIS® INTERNA HEX / Multi-Unit

NUCLE® OSS T6 BONELEVEL

## O series compatible with

OSTEOCARE® INTERNAL HEX

Osteo Ti Internal Hex

OSSTEM® GS(TS) / SS / SS OCTA Multi-Unit / US / Multi-Unit

OXY IMPLANT® K1 LINE / PSK

OXY K1 LINE / PSK Multi-Unit

OXY FIXO Multi-Unit

## P series compatible with

PHIBO® TSH

PLATON® PL-SUB

PALTOP® INTERNAL HEX / Multi-Unit

## Q series compatible with

Q IMPLANT® QZA

## R series compatible with

REX Piezol Implant TL 1.8 External / Multi-Unit

## S series compatible with

SGS DENTAL® INTERNAL HEXAGON / Multi-Unit

SHINHUNG®

SIC Invent® / Multi-Unit®

SIN® INTERNAL

SLH Implant DC Bone Level

SNUCONE Internal Submerged / Multi-Unit

SOUTHERN IMPLANTS® DEEP Conical / External / M-Series™ / IT Connection™ / Tri-Nex™ / Multi-Unit®

STRAUMANN® ANGLE Multi-Unit / BLX® / TLX® / Conical Bone Level® / synOcta® / Screw Retained Abutment® / Multi-Unit®

Sweden&Martina® CSR / Prama / Shelta / Premium Straight/SP Premium Kohno / Premium One / P.A.D.r Multi-Unit / P.A.D Multi-Unit

## T series compatible with

Tatum Implant

THOMMEN SPI®

TITANIUM FIX EXTERNA

TISMART2®

TRATE® Internal / Multi-Unit

Tree-Oss® Internal Hex

TRI® Narrow / Octa / Vent / Multi-Unit

## W series compatible with

WARANTEC® Oneplant

## Z series compatible with

ZIACOM® ZIACOM® / External

ZIMMER® Eztetic™ / Swiss Plus™ / Spline™ A / Spline™ B / Paragon™ / Tapered Screw-Vent® / Universal Abutment Multi-Unit / Tapered Screw-Vent® Multi-unit





THE LARGEST CONSUMABLE SYSTEM RANGE  
IN THE WORLD










**HOIL**<sup>3</sup> SHOP  
D E N T A L



[www.hoildentalshop.com](http://www.hoildentalshop.com)

Try Our  
Implant Prosthetic  
Components with Top Quality

Consumables

Consumables	Single	Pack
 <b>Digital Analog</b> 500+ free library setup Stone / 3D printed model	£10.00	Pack of 10 £79.95 (£7.99 each)
 <b>Screw - Straight</b> Raw Material: Ti 6Al-4V ELI (ASTM F 136) Same as the original screw design	£9.50	Pack of 10 £79.95 (£7.99 each)
 <b>Screw - Angled</b> Raw Material: Ti 6Al-4V ELI (ASTM F 136)	£12.00	Pack of 10 £99.95 (£9.99 each)
 <b>Intraoral / Model Scanbody</b> 900+ free library setup by remote support	£30.00	Pack of 5 £142.50 (£28.50 each)
 <b>Angled Ti-base</b> Raw Material: Ti 6Al-4V ELI (ASTM F 136) Inc. 2 Screws - Straight & Angled Screw  Up to 500+ Implant connections	£32.95	 <b>Straight Screw</b> Same as the original screw design  <b>Angled Screw</b>

VAT. Excluded.

iPen Lab Driver

Screw Driver Tips: **850+** Implant system available  
Analog Insertion Tips: **810+** Implant system available  
For the detailed guidance of iPen tips, please visit our online shop, [www.hoildentalshop.com](http://www.hoildentalshop.com).



Scan for the full 850+ system!



iPen Lab Driver

Price




iPen Set 1 - 10 Tips Package  
iPen Base + iPen Handle + 10 Screw Driver Tips

Special Price  
£278.00  
(Regular Price £330)



iPen Set 2 - 5 Tips Package  
iPen Base + iPen Handle +  
'5 Analog Insertion Tips' or 'Any 5 Screwdriver Tips'

Special Price  
£168.00  
(Regular Price £200)



iPen Base

£38.00



iPen Handle

£32.00



Screw Driver Tip (each)



£26.00



Analog Insertion Tip (each)

£26.00

For the detailed guidance of Clinical driver, please visit our online shop.

Type		Implant System	Recommended Torque
	Ball Hex	1. Some of angled screws of ARUM Ti-Base (screw head size should be less than 2.1mm) 2. ARUM Scanbody Screw	15-20 Ncm
	Ball Torx	1. Most of angled screws of ARUM Ti-Base (screw head size should be more than 2.1mm) 2. All implant system for disk angulated screw channels	15-20 Ncm

Clinical Ball Screwdriver

Compatible  
with  
Straumann®  
Ratchet

Order Code

Price



Driver Ball Hex - 15mm  
Ti-base Angled Screw/Scanbody



BH1215

£37.00



Driver Ball Hex - 22mm  
Ti-base Angled Screw / Scanbody



BH1222

£37.00



Driver Ball Torx - 15mm  
Ti-base Angled Screw/Angled Milling Screw



BT1715

£37.00



Driver Ball Torx - 22mm  
Ti-base Angled Screw/Angled Milling Screw



BT1722

£37.00



Driver Tip Ball Hex - 25mm  
Ti-base Angled Screw/Scanbody

BHD12H

£30.00











Driver Tip Ball Torx - 25mm  
Ti-base Angled Screw/Angled Milling Screw





BTD17H





£30.00







HOIL <sup>®</sup> SHOP Product List							
No	Compatible Product Group	Size	Angled Ti-Base		Model/ Intraoral Scanbody	Digital Analog	Screw
			H	NH			
1	ADIN <sup>®</sup> CLOSEFIT <sup>™</sup>	3.0	Ti_606	Ti_607_NH	B148	F606	S606
		3.5	Ti_610	Ti_611_NH	B149	F610	S610
		4.3/5.0	Ti_614	Ti_614_NH	B150	F616	S614
2	ADIN <sup>®</sup> TOUAREG <sup>™</sup> S&OS&SWELL	3.5/3.75/4.2/5.0/6.0	Ti_168	Ti_598_NH	B081	F168	S168
3	ADIN <sup>®</sup> Multi-Unit	-	-	-	M218_00-ADIN	MF200-ADIN	MS244
4	ADIN <sup>®</sup> UnitFit <sup>™</sup>	-	-	-	B129-ADIN	-	Contact us
5	ALPHA-BIOTEC <sup>®</sup> CONICAL	3.2/3.25/3.5	-	-	B151-ALP	-	Contact us
		3.75/4.2/5.0	-	-	B016_00-ALP	-	Contact us
6	ALPHA-BIO TEC <sup>®</sup> Internal Hex	3.75/4.2/5.0/6.0	Ti_169	Ti_600_NH	B082	F169	S169
7	ALPHA-BIO TEC <sup>®</sup> Multi-Unit	TCT-N	-	-	M237-ALP	-	Contact us
8	ANTHOGYR <sup>®</sup> AnthoFit <sup>®</sup> OI	3.5	-	-	B114	-	S229
		4.0	-	-	B114	-	S229
		5.0	-	-	B114	-	S229
9	ANTHOGYR <sup>®</sup> Axiom <sup>®</sup> BL	-	Ti_137	Ti_365_NH	B080	F137	S137
10	ANTHOGYR <sup>®</sup> Axiom <sup>®</sup> TL	4.0	-	-	-	-	Contact us
		4.8	-	-	-	-	Contact us
11	ANTHOGYR <sup>®</sup> (Axiom <sup>®</sup> ) Multi-Unit <sup>®</sup>	MU-40	-	-	M223	-	MS222
		MU-48	-	-	M222	-	MS222
12	ASTRA TECH <sup>™</sup> EV	3.0	Ti_269	Ti_508_NH	B128	F269	S269
		3.6	Ti_270	Ti_510_NH	B129	F270	S270
		4.2	Ti_271	Ti_512_NH	B130	F271	S271
		4.8	Ti_272	Ti_514_NH	B131	F272	S272
		5.4	Ti_273	Ti_516_NH	B132	F272	S272
13	ASTRA TECH <sup>™</sup> EV <sup>™</sup> PROFILE	4.2	-	-	B138	-	S271
		4.8	-	-	B141	-	S272
14	ASTRA TECH <sup>™</sup> OsseoSpeed <sup>™</sup> TX	YELLOW 3.0	Ti_046	Ti_335_NH	B023	F046	S046
		AQUA 3.5/4.0	Ti_047	Ti_048_NH	B024	F047	S047
		LILAC 4.5/5.0	Ti_049	Ti_050_NH	B025	F049	S049
15	ASTRA TECH <sup>™</sup> Multi-Unit <sup>®</sup>	UNI20	-	-	M209	-	MS209
		EV UNI30	-	-	M242	-	MS672
		UNI45	-	-	M210	-	MS209
		Angled(TCA)	-	-	M211	-	MS209
16	ASTRA TECH <sup>™</sup> EV <sup>™</sup> Multi-Base	-	-	-	M200_00-EVM	-	Contact us
17	AVINENT <sup>®</sup> OCEAN	-	-	-	B093-AV	F189-AV	Contact us
18	B&B DENTAL Conexa	3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0	-	-	B115-BB	-	Contact us
19	B&B <sup>®</sup> Multi-Unit <sup>®</sup>	-	-	-	M221-BB	-	Contact us
20	BEGO <sup>®</sup> Internal	3.0	-	-	-	-	Contact us
		3.25/3.75	Ti_053	Ti_518_NH	B026	F053	S053
		4.1	Ti_054	Ti_520_NH	B027	F054	S053
		4.5	Ti_055	Ti_522_NH	B027	F055	S053
		5.5	Ti_056	Ti_524_NH	B027	F056	S053
21	BEGO <sup>®</sup> Multi-Unit <sup>®</sup>	-	-	-	M208	-	MS208





www.hoildentalshop.com							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
22	BICON <sup>®</sup>	2.0	Ti_057	-	B028	F057	-
		2.5	Ti_058	-	B029	F058	-
		3.0	Ti_059	-	B030	F059	-
23	BIO3 IMPLANT Conical Connection	Standard 3.3/3.8	-	-	B015-B3CC	-	Contact us
		Wide 4.2/5.0	-	-	B016_00-B3CC	-	Contact us
24	BIOCOMP <sup>®</sup> Internal Hex	4.0	-	-	-	-	-
25	BIOCOMP <sup>®</sup> Multi-Unit	4.0	-	-	M251-BIOC	-	-
26	BIODENTA <sup>®</sup> Bone Level	B1 3.5	-	-	B178	F663	S683
		B2 4.1/4.8/6.0	-	-	B179	F664	S684
27	BIOHORIZONS <sup>®</sup> External <sup>®</sup>	3.5(NP)	Ti_192	Ti_193_NH	B017_02	F192	S192
		4.0(RP)	Ti_194	Ti_195_NH	B018_05	F194	S194
		5.0(WP)	Ti_196	Ti_197_NH	B098	F196	S194
		6.0	Ti_198	Ti_199_NH	B098	F198	S194
28	BIOHORIZONS <sup>®</sup> Internal <sup>®</sup>	3.0	Ti_291	Ti_292_NH	B140	F291	S291
		3.5(NP)	Ti_189	Ti_526_NH	B093	F189	S291
		4.5(RP)	Ti_190	Ti_528_NH	B094	F190	S291
		5.7(WP)	Ti_191	Ti_530_NH	B095	F191	S291
29	BIOHORIZONS <sup>®</sup> Tissue level	3.5	-	-	B093-BIOH	-	Contact us
		4.5	-	-	B094-BIOH	-	Contact us
		5.7	-	-	B095-BIOH	-	Contact us
30	BIOHORIZONS <sup>®</sup> Multi-Unit <sup>®</sup>	4.8	-	-	M200_00-BIOH	MF200-BIOH	MS245
31	BIOMET 3i <sup>®</sup> Certain <sup>®</sup>	3.4	Ti_065	Ti_534_NH	B033	F065	S065
		4.1	Ti_066	Ti_536_NH	B034	F066	S065
		5.0	Ti_067	Ti_538_NH	B035	F067	S065
		6.0	Ti_068	Ti_540_NH	B036	F068	S065
32	BIOMET 3i <sup>®</sup> External <sup>®</sup>	MINI	Ti_060	Ti_060_NH	B031_00	F060	S060
		RP	Ti_061	Ti_062_NH	B018_01	F061	S061
		WP	Ti_063	Ti_064_NH	B018_01	F063	S061
		XWP	-	-	B018_01	-	S061
33	BIOMET 3i <sup>®</sup> Multi-Unit <sup>®</sup>	-	-	-	M218_00	-	MS218
34	BIONER IMPLANT	3.5	-	-	B054-BIO	-	Contact us
35	BIOTEC Internal Hex	-	-	-	B093-BTC	-	Contact us
36	BIOTEC Multi-Unit	-	-	-	M228-BTC	-	Contact us
37	BIOTECH <sup>®</sup> BIS EXTERNAL	3.6	-	-	-	-	Contact us
38	BIOTECH <sup>®</sup> BIS Octagone Interne	3.9	-	-	B114-BT	-	Contact us
		4.4	-	-	B114-BT	-	Contact us
		5.4	-	-	B114-BT	-	Contact us
39	BIOTECH <sup>®</sup> Kontact	3.0	-	-	B173	F661	S679
		3.6/4.2/4.8/5.4	Ti_201	Ti_202_NH	B100	F201	S201
40	BIOTECH <sup>®</sup> Multi-Unit	Narrow/4.0	-	-	M218_00-BT	-	Contact us
		4.9	-	-	M218_00-BT	-	Contact us
41	BrainBase Corporation <sup>®</sup> MYTIS <sup>®</sup> Arrow Implant	B Type	-	-	B109	-	S215
		C Type	-	-	B109	-	S215
		E Type	-	-	B109	-	S215
		EW Type	-	-	B112	-	S215





HOIL <sup>®</sup> SHOP Product List							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
42	BREDENT Medical Copa SKY®	4.0/5.0/6.0	-	-	-	-	Contact us
43	BREDENT Medical SKY®	3.5	Ti_290	Ti_542_NH	B139	F290	S290
		4.0	Ti_284	Ti_544_NH	B135	F284	S290
		MINI2	-	Ti_546_NH	B136	F285	S285
44	BREDENT Medical SKY® Multi-Unit®	UNI.CON	-	-	M228	-	MS229
		FAST&FIXED	-	-	M229	-	MS229
45	BTI® EXTERNA®	3.5	-	-	B031_00-BTI	-	Contact us
		4.1	-	-	B018_01-BTI	F061-BTI	Contact us
		5.5	-	-	B019_01-BTI	-	Contact us
46	BTI® INTERNA®	3.5 MINI	Ti_646	-	B154	F646	S646
		4.1 REGULAR	Ti_647	-	B155	F647	S646
		5.5 WIDE	Ti_651	-	B155	F651	S646
47	BTI® Multi-Unit®	3.5	-	-	M209-BTI	-	MS241
		4.1	-	-	M241	MF241	MS241
		5.5	-	-	M241	MF241	MS241
		5.0 ANG CONVERTER	-	-	M221-BT	-	Contact us
48	BTK® Multi-Unit	-	-	-	M200_00-BTK	-	Contact us
49	CAMLOG®	3.3	Ti_069	Ti_548_NH	B037	F069	S069
		3.8	Ti_070	Ti_550_NH	B038	F070	S069
		4.3	Ti_071	Ti_552_NH	B039	F071	S069
		5.0	Ti_072	Ti_554_NH	B040	F072	S072
50	CAMLOG® iSy®	-	-	-	B015-CAM	-	Contact us
51	CAMLOG® Multi-Unit	3.3/3.8/4.3	-	-	M251	MF242	MS251
		5.0/6.0	-	-	M252	MF243	MS252
52	CONELOG®	3.3	Ti_295	Ti_556_NH	B144	F295	S295
		3.8	Ti_296	Ti_558_38_NH	B145	F296_38	S295
		4.3	Ti_296	-	B145	F296_43	S295
		5.0	Ti_297	Ti_560_NH	B146	F297	S297
53	CONELOG® Multi-Unit	3.3/3.8/4.3	-	-	M251	MF242	MS251
		5.0/6.0	-	-	M252	MF243	MS252
54	CHAMPION® Multi-Unit	-	-	-	-	-	-
55	CLC conic™	3.5/4.0/4.5/5.0/6.0	-	-	B221	-	S622
56	CONEXAO HI (Hexagono Interno)	4.1	-	-	B043-CON	-	Contact us
57	CORE VENT® Internal	3.25 Yellow	-	-	-	-	Contact us
58	CORTEX®	3.3/3.8/4.2/5.0/6.0	Ti_265	-	B125	F265	S265
59	CORTEX® Multi-Unit	-	-	-	M232	-	MS232
60	COWELLMEDI® INNO Internal	4.0/4.5/5.0/6.0	Ti_001	Ti_323_NH	B001_00	F001	S001
61	C-TECH® Bone level	3.8/4.3/5.1	Ti_230	Ti_639_NH	B115	F230	Contact us
62	C-TECH® Esthetic Line	3.8/4.3/5.1	Ti_230	Ti_639_NH	B115	F230	S230
63	C-TECH® Multi-Unit	-	-	-	M221-CTE	-	Contact us
64	DEEP®	3.8 H	Ti_279	Ti_288_NH	B127	F279	S279
		4.5 H	Ti_280	Ti_289_NH	B127	F280	S279
65	DENTAURUM® TIO	3.3	-	-	B167	-	-
		3.7/4.2	-	-	B162	-	-
		4.8/5.5	-	-	B163	-	-





www.hoildentalshop.com							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
66	DENTIS® i-clean	4.8	Ti_002	Ti_449_NH	B002	F002	S002
		6.5	Ti_237	Ti_441_NH	B061_00	F237	S237
67	DENTIS® Save	4.8	-	-	B002-DEN	-	Contact us
		6.5	-	-	B061_00-DEN	-	Contact us
68	DENTIS® s-clean	Narrow	Ti_181	Ti_325_NH	B085	F181	S181
		Regular/Wide	Ti_003	Ti_003_NH	B001_01	F003	S003
69	DENTIS® Multi-Unit	-	-	-	M234	-	MS234
70	DENTIUM® NR LINE	NR	Ti_004	Ti_368_NH	B003	F004	S004
		NR36	Ti_236	Ti_506_NH	B117	F236	S004
71	DENTIUM® Simple Line	4.8	Ti_128	Ti_451_NH	B156	F128	S128
		6.5	Ti_129	Ti_453_NH	B157	F129	S128
72	DENTIUM® Super Line	3.6/4.0/4.5/5.0/6.0/7.0	Ti_005	Ti_007_NH	B001_02	F005	S005
73	DENTIUM® Multi-Unit	CONVERTIBLE 4.5	MTi_212	-	M212	MF212	MS212
		CONVERTIBLE 5.5	-	-	M239	MF239	MS212
74	DENTSPLY® Ankylos®	PILAR REGULAR 3.5/4.5/5.5/7.0	Ti_136	Ti_223_NH	B079	F136	S136
75	DENTSPLY® Ankylos® Balance base	Balance Base Narrow	MTi_235	-	M235	MF235	MS235
		Balance Base	MTi_236	-	M236	MF236	MS236
76	DENTSPLY® Xive®	3.0	Ti_073	Ti_562_NH	B041	F073	S073
		3.4	Ti_074	Ti_564_NH	B042	F074	S074
		3.8	Ti_075	Ti_566_NH	B043	F075	S074
		4.5	Ti_076	Ti_568_NH	B044	F076	S074
		5.5	Ti_077	Ti_570_NH	B045	F077	S074
77	DENTSPLY® Xive® Multi-Unit	3.4	-	-	M221-XIVE	-	Contact us
		3.8	-	-	M221-XIVE	-	Contact us
		4.5	-	-	M221-XIVE	-	Contact us
		5.5	-	-	M221-XIVE	-	Contact us
78	DIO® SM	Mini	Ti_009	Ti_357_NH	B004	F009	S009
		Regular/Wide/Extra Wide	Ti_011	Ti_359_NH	B005	F011	S011
79	DIO® SM	AMI 48	Ti_014	Ti_464_NH	B006	F014	S014
80	DIO® UF	Narrow	Ti_210	Ti_310_NH	B084_00	F210	S210
		Regular/Wide	Ti_013	Ti_278_NH	B001_03	F013	S013
		Non Submerged	-	-	B171	-	S676
81	DIO® (UF) Multi-Unit	MU-48	-	-	M227	MF227	MS227
82	DYNA® PUSHIN Octalock®	3.6/4.0/5.0	Ti_275	-	B126	F275	S275
83	EBI®	Octa ST	-	-	B104	-	S206
		Octa C1	-	-	B104	-	S206
		Octa C2	-	-	B104	-	S206
		Octa C3	-	-	B104	-	S206
		TL48	-	-	B060_05	-	S206
84	EDISON MEDICAL® Clicq Multi-Unit	-	-	-	M200_01-EDI	-	Contact us
85	EUROTEKNIKA™ Aesthetica	4.8	-	-	B060_00-EURO	F101-EURO	Contact us







HOIL <sup>®</sup> SHOP Product List							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
86	EUROTEKNIKA™ Naturactis / Naturall+ / Natea+	3.0	-	-	B084_00-EURO	-	Contact us
		3.5/4.0/4.5/5.0	-	-	B024-EURO	-	Contact us
87	EUROTEKNIKA™ Multi-Unit	4.8	-	-	M218_00-EURO	MF200-EURO	Contact us
88	GC IMPLANT® Aadva	Narrow	-	-	-	-	Contact us
		Regular	-	-	B011_00-GC	-	Contact us
89	Global D® Tekka® In-kone	3.5/4.0/4.5/5.5	Ti_203	Ti_308_NH	B101	F203	S203
90	Global D® In-Kone® Multi-Unit	STRAIGHT	-	-	M215	-	MS215
		ANGLED	-	-	M216	-	MS216
91	Global D® EVL® Multi-Unit	ST/WD	-	-	M227-GD	-	Contact us
92	GP Implant Multi-Unit	-	-	-	M200_00-GP	-	Contact us
93	HIOSSSEN® ET	Mini 3.5	Ti_029	Ti_031_NH	B015	F029	S029
		Standard 4.0/4.5/5.0	Ti_032	Ti_033_NH	B016_00	F032	S032
94	HIOSSSEN® SS	4.8(Regular)	Ti_043	Ti_043_NH	B020	F043	S043
		6.0(Wide)	Ti_044	Ti_044_NH	B021	F044	S043
95	HIOSSSEN® Multi-Unit	Mini	-	-	M213	-	MS213
		Regular	-	-	M213	-	MS213
96	HumanTech RATIO®	MINI	-	-	B222	-	S628
		STANDARD	-	-	B223	-	S628
		LARGE	-	-	B224	-	S628
97	IHDEDENTAL Hexacone®	-	-	-	B093-IHD	-	Contact us
98	IMPLANT DIRECT® InterActive™	3.0	-	-	B054-ID	F090-ID	Contact us
		3.4	-	-	B055-ID	F091-ID	Contact us
99	IMPLANT DIRECT® Legacy™	3.0	Ti_224	Ti_349_NH	B113	F224	S224
		3.5	Ti_212	Ti_572_NH	B062_01	F212	S212
		4.5	Ti_213	Ti_574_NH	B063_01	F213	S212
		5.7	Ti_214	Ti_576_NH	B064_01	F214	S212
100	IMPLANT DIRECT® Swish Active™	3.0/3.3	-	-	B054-ID	F090-ID	Contact us
		3.4/4.1/4.8	-	-	B055-ID	F091-ID	Contact us
101	IMPLANT DIRECT® Swish Plus™	4.8	-	-	B060_00-ID	F101-ID	Contact us
		6.5	-	-	B061_01-ID	F103-ID	Contact us
102	IMPLANT DIRECT® Tri-Lobe™	3.5	-	-	B050-ID	F082-ID	Contact us
		4.3	-	-	B051-ID	F084-ID	Contact us
		5.0	-	-	B052-ID	F086-ID	Contact us
103	IMPLANT DIRECT® Multi-Unit	-	-	-	M214	-	MS214
104	IMPLANT SWISS® Bone Level	3.3	-	-	B024-IPS	-	Contact us
		3.7/4.3/4.8/5.5	-	-	B016_00-IPS	-	Contact us
105	IMPLANT SWISS Multi-Unit	4.8	-	-	M200_00-IPS	-	Contact us
106	INTRA LOCK GOLD&BLUE	STANDARD	-	-	B165	-	S672
		WIDE	-	-	B164	-	S672
107	IPD Multi-Unit	-	-	-	M200_00-IPD	-	Contact us
108	iRES® Internal Hex 2.5mm	3.7/4.1/4.7/5.2	-	-	B062_00-IRS	-	Contact us
109	iRES® Internal Octa Tissue	3.7/4.1/4.7	-	-	B060_00-IRS	-	Contact us
110	iRES® Multi-Unit	-	-	-	M222-IRS	-	Contact us
111	J DENTAL CARE® JD Evolution®	S	-	-	B090-JD	-	Contact us
		3.7/4.3/5.0/6.0	-	-	B026-JD	-	Contact us





www.hoildentalshop.com							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
112	J DENTAL CARE® JD Evolution® Plus+	3.7/4.3/5.0/6.0	-	-	B093-JD	F189-JD	Contact us
113	J DENTAL CARE® JD Icon®	3.9/4.3/5.0	-	-	B055-JD	F091-JD	Contact us
114	J DENTAL CARE® JD Icon® Plus+	3.7/4.3/5.0	-	-	B055-JD	F091-JD	Contact us
115	J DENTAL CARE® JD OCTA®	4.8	-	-	B060_00-JD	F101-JD	Contact us
		6.5	-	-	B061_01-JD	F103-JD	Contact us
116	J DENTAL CARE® Multi-Unit	-	-	-	M200_00-JD	MF200-JD	MS247
117	KENTEC®	SB1	Ti_241	-	B122	F241	S241
		SB2	Ti_242	-	B123	F242	S241
		SB3	Ti_243	-	B124	F243	S243
118	KEYSTONE® Prima Connex®	3.5	Ti_244	-	B116	F244	S244
		4.1	Ti_245	-	B119	F245	S244
		5.0	Ti_246	-	B120	F246	S244
119	Klockner® SK2	-	-	-	-	-	Contact us
120	KYOCERA® FINESIA® Bone-level	3.7/4.2	-	-	B168	-	S673
		4.7/5.2	-	-	B169	-	S674
121	KYOCERA® Poiox	3.4	Ti_211	Ti_652_NH	B105	-	S211
		3.7	Ti_078	Ti_653_NH	B046	F078	S078
		4.2	Ti_079	Ti_654_NH	B047	F079	S079
		4.7	Ti_080	Ti_655_NH	B048	-	S080
		5.2	Ti_081	Ti_656_NH	B049	-	S080
122	LASAK Bioniq®	S3.5/S4.0/S5.0/T4.0/T5.0	Ti_204	Ti_578_NH	B102	F204	S204
123	LASAK ® Multi-Unit	-	-	-	M200_00-LAS	-	-
124	LEADER Italia® Mores Connection	-	-	-	B024-LDI	-	Contact us
125	MEGAGEN® AnyONE	3.5/4.0/4.5/5.0/5.5/6.0/7.0	Ti_023	Ti_024_NH	B012	F023	S023
126	MEGAGEN® AnyONE MINI	3.0	Ti_170	Ti_282_NH	B084_01	F170	S170
127	MEGAGEN® Anyridge	Small 3.5/R/W/SW	Ti_021	Ti_131_NH	B011_00	F021	S021
		EW 3.3/4.0/4.8	-	-	B011_01	F021	S641
128	MEGAGEN® Blue Diamond	NC	-	-	-	-	Contact us
		RC	-	-	-	-	Contact us
129	MEGAGEN® EZ Plus	Mini	Ti_015	Ti_130_NH	B007	F015	S015
		Regular/ Wide	Ti_016	Ti_017_NH	B008	F016	S016
130	MEGAGEN® Exfell	4.8	Ti_267	Ti_445_NH	B060_03	F267	S267
131	MEGAGEN® MegaFix	-	Ti_025	Ti_328_NH	B013	F025	S025
132	MEGAGEN® RESCUE	External 5.0	Ti_020_50	Ti_401_NH_50	B010	F020_50	S020
		External 6.0	Ti_020_60	Ti_401_NH_60	B010	F020_60	S020
		Internal	Ti_019	Ti_248_NH	B009	F019	S019
133	MEGAGEN® Multi-Unit	S-Type 4.8	-	-	M221-MG	-	Contact us
		OCTA 3.8	-	-	-	-	Contact us
		OCTA 4.8	-	-	M202-MG	-	Contact us
		N-Type 4.8	-	-	M200_00-MEG	MF200-MEG	MS246
		RC-MUA	-	-	M218_00-MEG	-	Contact us
134	MIS® C1	NARROW	Ti_183	Ti_633_NH	B087	F183	S183
		STANDARD	Ti_184	Ti_381_NH	B088	F184	S184





HOIL <sup>®</sup> SHOP Product List							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
134	MIS® C1	WIDE	Ti_185	Ti_384_NH	B089	F185	S184
135	MIS® Connect	-	-	-	B090-MIS	-	Contact us
136	MIS® Internal Hexagon	NARROW	Ti_186	Ti_636_NH	B090	F186	S186
		STANDARD	Ti_187	Ti_637_NH	B091	F187	S187
		WIDE	Ti_188	Ti_638_NH	B092	F188	S187
137	MIS® V3	NARROW	-	-	-	-	Contact us
		STANDARD	Ti_184-MIS	Ti_381_NH-MIS	B088-MIS	F184-MIS	Contact us
138	MIS® Multi-Unit	STANDARD	-	-	M217	-	MS217
139	MIS® Multi-Unit (NEW)	Multi (New)	-	-	M218_00-MIS	-	Contact us
140	MSDI Alef-Spiral Internal Hex	3.5/3.75/4.2/5.0/6.0	-	-	B093-AL	-	Contact us
141	Medentika® Multi-Unit	-	-	-	M200_00-MDT	MF200-MDT	Contact us
142	Medentis Medical® ICX®	3.75/4.1/4.8	Ti_205	Ti_353_NH	B103	F205	S205
143	Medentis Medical® ICX® ACTIVE Master®	3.3	-	-	B015-MMIC	-	Contact us
144	Medentis Medical® ICX® Multi-Unit	ICX-INS	-	-	M200_01	-	MS233
145	NEOBIOTECH® IS ACTIVE SCRP	-	Ti_268	-	B001_04	F026	S026
146	NEOBIOTECH® IS SYSTEM	-	Ti_026	Ti_027_NH	B001_04	F026	S026
147	NEOBIOTECH® IT SYSTEM	-	Ti_226	-	B060_04	F226	S226
148	NEODENT® CM	3.5/3.75/4.0/4.3/5.0	Ti_643	-	B151	F643	S643
149	NEODENT® GM	3.5/3.75/4.0/4.3/5.0/6.0	-	-	B152	F644	S644
150	NEODENT® External	3.3	-	-	B017_02-NEO	-	Contact us
		4.1	-	-	B018_04-NEO	-	Contact us
		4.3	-	-	B018_04-NEO	-	Contact us
		5.0	-	-	B019_01-NEO	-	Contact us
151	NEODENT® CM Conical	Micro	-	-	M217-NED	-	MS243
		Mini	-	-	M218_00-NED	MF200-NEOC	MS243
152	NEODENT® GM Conical	Micro	-	-	M217-GM	-	Contact us
		Mini	-	-	M218_00-GM	MF200-NEOG	Contact us
153	NEOSS Access®	-	-	-	-	-	Contact us
154	NEOSS Pro Active®	3.25	-	Ti_618_NH	B219	F474	S474
		3.5	-	Ti_620_NH	B220	F476	S476
155	NOBEL BIO CARE® Active™	3.0	Ti_281	Ti_580_NH	B134	F281	S281
		3.5(NP)	Ti_090	Ti_341_NH	B054	F090	S090
		4.3/5.0(RP)	Ti_091	Ti_343_NH	B055	F091	S091
		WP	Ti_266	Ti_582_NH	B133	F266	S091
156	NOBEL BIO CARE® Branemark®	3.5(NP)	Ti_118	Ti_119_NH	B017_01	F118	S118
		4.0(RP)	Ti_120	Ti_121_NH	B018_04	F120	S120
		5.0(WP)	Ti_122	Ti_123_NH	B019_01	F122	S122
157	NOBEL BIO CARE® N1™	3.5(NP)	-	-	-	-	-
		4.0/4.8(RP)	-	-	-	-	-
158	NOBEL BIO CARE® Replace Select™	3.5(NP)	Ti_082	Ti_584_NH	B050	F082	S082
		4.3(RP)	Ti_084	Ti_586_NH	B051	F084	S084

www.hoildentalshop.com							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
158	NOBEL BIO CARE® Replace Select™	5.0(WP)	Ti_086	Ti_588_NH	B052	F086	S084
		6.0(SW)	Ti_088	Ti_590_NH	B053	F088	S084
159	NOBEL BIO CARE® MULTI AESTHETIC	-	-	-	-	-	Contact us
160	NOBEL BIO CARE® Multi-Unit™	4.8	MTi_200	-	M200_00	MF200	MS200_00
		6.0	MTi_201	-	M201	MF201	MS201
		Abutment Plus	MTi_680	-	M218_01	MF680	MS200_01
161	NORIS® INTERNAL HEX	3.75/4.2/5.0/6.0	-	-	B062_00-NOR	-	Contact us
162	NORIS® Multi-Unit	-	-	-	M200_00-NOR	-	Contact us
163	NUCLE® OSS T6 BONELEVEL	Nr/Sd/Wd	-	-	B159	-	S669
164	OSSTEM® GS (TS)	Mini	Ti_029	Ti_031_NH	B015	F029	S029
		Regular/Ultra-Wide	Ti_032	Ti_033_NH	B016_00	F032	S032
165	OSSTEM® SS	4.8(Regular)	Ti_043	Ti_043_NH	B020	F043	S043
		6.0(Wide)	Ti_044	Ti_044_NH	B021	F044	S043
166	OSSTEM® SS OCTA Multi-Unit	4.8(Regular)	-	-	M202-OSS-MT	-	Contact us
167	OSSTEM® US	3.5(Mini)	Ti_036	Ti_037_NH	B017_00	F036	S036
		4.1(Regular)	Ti_038	Ti_039_NH	B018_00	F038	S038
		5.1(Wide)	Ti_041	Ti_042_NH	B019_00	F041	S041
168	OSSTEM® Multi-Unit	Mini	-	-	M213	-	MS213
		Regular	-	-	M213	-	MS213
		TS-CON-40	-	-	M224	-	MS224
		TS-CON-50	-	-	M225	-	MS225
		TS-CON-60	-	-	M226	-	MS225
169	OSTEOCARE® INTERNAL HEX	3.75	-	-	-	-	Contact us
		4.5	-	-	-	-	Contact us
170	OXY Implant® K1 Line	Mini	-	-	-	-	Contact us
		Regular	-	-	B011_00-OXY	-	Contact us
171	OXY Implant® PSK	Mini	-	-	-	-	Contact us
		Standard	-	-	B011_00-OXY	-	Contact us
172	OXY Implant® K1 Line/PSK Multi-Unit	-	-	-	M221-OXY	-	Contact us
173	OXY Implant® FIXO Multi-Unit	Mini	-	-	M207-OXY	-	Contact us
		Regular	-	-	M221-OXYF	-	Contact us
174	Osteo Ti Internal Hex	-	-	-	-	-	Contact us
175	PALTOP® INTERNAL HEX	3.25	-	-	B090-PAL	-	Contact us
		3.75/4.2/5.0	-	-	B082-PAL	-	Contact us
176	PALTOP® Multi-Unit	NP	-	-	-	-	Contact us
		RP	-	-	-	-	Contact us
177	PHIBO® TSH	4.0	-	-	B018_01-PH	-	Contact us
178	PLATON® PL-SUB	3.3/3.8/4.7	-	-	B086	-	S182
179	Q IMPLANT® QZA	3.5	-	-	B091-QI	-	Contact us
180	REX Piezo Implant TL 1.8 External	-	-	-	B018_01-REX	-	Contact us
181	REX IMPLANTS® Multi-Unit	-	-	-	M237-REX	-	Contact us
182	SGS Dental® Internal Hexagon	3.5/3.75/4.2/4.5 /5.0/6.0	-	-	B093-SGS	-	Contact us
183	SGS DENTAL® Multi-Unit	-	-	-	M218_00-SGS	-	Contact us



HOIL <sup>®</sup> SHOP Product List							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
184	SHINHUNG®	M	-	-	B022	-	S045
		LUNA NP	-	-	B015-SHL	-	Contact us
		LUNA RP	-	-	B016_01	F274	S274
185	SIC Invent®	3.3	Ti_092	-	B056	F092	S092
		4.2	Ti_094	-	B057	F094	S092
186	SIC Multi-Unit	-	-	-	M221	-	MS221
187	SIN IMPLANT® Internal	-	-	-	B151-SN	-	Contact us
188	SLH Implant DC Bone Level	-	-	-	-	-	Contact us
189	SNUCONE® Internal Submerged	3.5	-	-	B015-SNU	-	Contact us
		4.0/4.3/4.8/5.3/5.8/6.3/6.8	-	-	B001_02-SNU	-	Contact us
190	SNUCONE® Multi-Unit	4.5	-	-	M200_00-SNC	-	Contact us
191	SOUTHERN IMPLANTS® DEEP Conical	3.0	Ti_046-SIDC	Ti_335_NH-SIDC	B023-SIDC	F046-SIDC	Contact us
		3.5/4.0	Ti_047-SIDC	Ti_048_NH-SIDC	B024-SIDC	F047-SIDC	Contact us
		4.5/5.0	Ti_049-SIDC	Ti_050_NH-SIDC	B025-SIDC	F049-SIDC	Contact us
192	SOUTHERN IMPLANTS® External	3.25	Ti_626	Ti_627_NH	B031_01	F626	S631
		4.0	Contact us	Contact us	B018_02	F061-SIDC	S631
		5.0	Contact us	Contact us	B018_02	F063-SIDC	S631
		6.0	Ti_631	Ti_632_NH	B018_02	F631	S631
193	SOUTHERN IMPLANTS® IT Connection™	4.8	-	-	B060_00-SIDC	-	Contact us
		6.5	-	-	B061_01-SIDC	-	Contact us
194	SOUTHERN IMPLANTS® M-Series™	-	-	-	B166	F658	S677
195	SOUTHERN IMPLANTS® Tri-Nex™	3.5	-	-	B050-SIDC	F082-SIDC	Contact us
		4.3	-	-	B051-SIDC	F084-SIDC	Contact us
		5.0	-	-	B052-SIDC	F086-SIDC	Contact us
		6.0	-	-	B053-SIDC	F088-SIDC	Contact us
196	SOUTHERN IMPLANTS® Multi-Unit	4.8	-	-	M200_00-SIDC	MF200-SIDC	MS242
		6.0	-	-	M201-SIDC	-	MS242
197	STRAUMANN® BLX®	3.5/3.75/4.0/4.5/5.0/5.5/6.0(RB)	-	-	-	-	Contact us
		5.0/5.5/6.5(WB)	-	-	-	-	Contact us
198	STRAUMANN® TLX	3.75/4.5(NT)	-	-	-	-	Contact us
		3.75/4.5(RT)	-	-	-	-	Contact us
		5.5/6.5(WT)	-	-	-	-	Contact us
199	STRAUMANN® Conical Bone Level®	2.9(SC)	-	-	B172	F660	S678
		3.3(NC)	Ti_498	Ti_501_NH	B058_01	F096_01	S498
		4.1/4.8(RC)	Ti_499	Ti_504_NH	B059_01	F098_01	S499
200	STRAUMANN® synOcta®	3.5(NNC)	Ti_602	Ti_603_NH	B147	F602	S602
		3.5(NN)	Ti_239	Ti_428_NH	B121	F239	S239
		4.8(RN)	Ti_101	Ti_455_NH	B060_00	F101	S101
		6.5(WN)	Ti_103	Ti_457_NH	B061_01	F103	S101

www.hoildentalshop.com							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
201	STRAUMANN® ANGLE Multi-Unit	NC 25	-	-	M207	-	MS207
		RC 25	-	-	M207	-	MS207
202	STRAUMANN® Screw Retained Abutment Multi-Unit	3.5	MTi_237	-	M237	MF237	MS237
		4.6	MTi_238	-	M238	MF238	MS237
203	STRAUMANN® Multi-Unit	(synOcta) RN	-	-	M202	-	MS202
		(synOcta) WN	-	-	M203	-	MS202
		(Bone Level) 3.5(NC)	-	-	M204	MF204	MS204
		(Bone Level) 4.5(NC)	-	-	M205	MF205_NC	MS204
		(Bone Level) 4.5(RC)	-	-	M205	MF205_RC	MS204
		(Bone Level) 6.5(RC)	-	-	M206	MF206	MS204
204	Sweden&Martina® CSR	Small 3.0/3.5	-	-	B151-SW	-	Contact us
		Regular 3.8/4.2/5.0	-	-	B103-SW	-	S668
205	Sweden&Martina® Prama	3.3/3.8/4.25/5.0	Ti_670	-	B160	F652	S668
206	Sweden&Martina® Premium Kohno	3.3	Ti_670	-	B160	F652	S668
		3.8	Ti_671	-	B158	F653	S668
		4.25	Ti_672	-	B161	F654	S670
		5.0	-	-	B161	-	S670
		6.0	-	-	B161	-	S670
207	Sweden&Martina® Premium Straight/SP	3.3	Ti_670	-	B160	F652	S668
		3.8	Ti_671	-	B158	F653	S668
		4.25	Ti_672	-	B161	F654	S670
		5.0	-	-	B161	-	S670
208	Sweden&Martina® Premium one	3.3	Ti_670	-	B160	F652	S668
		3.8	Ti_671	-	B158	F653	S668
		4.25	-	-	B158	-	S668
		5.0	-	-	B158	-	S668
209	Sweden&Martina® Shelta	3.8	Ti_671	-	B158	F653	S668
		4.25	-	-	B158	-	S668
		5.0	-	-	B158	-	S668
		6.0	-	-	B158	-	S668
210	Sweden&Martina® PAD Multi-Unit	-	-	-	M221-SM	-	Contact us
211	Sweden&Martina® PAD.r Multi-Unit	-	-	-	M200_00-SM	MF200-SM	Contact us
212	THOMMEN SPI®	3.5	Ti_113	-	B069	F113	S113
		4.0	Ti_114	-	B070	F114	S114
		4.5	Ti_115	-	B071	F115	S114
		5.0	Ti_116	-	B071	F116	S114
		6.0	Ti_117	-	B071	F117	S114
213	TISMART2 Internal	3.5/3.8/4.5/5.0	-	-	B024-TIS	-	Contact us
214	TITANIUM FIX EXTERNA	3.5/4.1/5.1	-	-	B017_00-TIFIX	-	Contact us
215	TRATE® Internal	-	-	-	B151-TRT	-	Contact us
216	TRATE® Multi-unit	Small	-	-	-	-	Contact us
217	TRI® Narrow	3.3	-	-	-	-	S675
218	TRI® Octa	4.8	-	-	B060_00-TRI	F101-TRI	Contact us
219	TRI® Vent	3.5	-	-	B062_00-TRI	-	Contact us

<div> <div> <div>HOIL<sup>®</sup></div> <div>SHOP</div> </div> <div>Product List</div> </div>							
No	Compatible Product Group	Size	Angled Ti-Base		Intraoral/ Model Scanbody	Digital Analog	Screw
			H	NH			
220	TRI® Multi-Unit	-	-	-	M251-TRI	-	Contact us
221	Tree-Oss® Internal Hex	-	-	-	B093-TRO	-	Contact us
222	WARANTEC Oneplant	Tap4.3/5.3/ St3.6/4.1/5.1	Ti_028	Ti_433_NH	B014	F028	S028
223	ZIACOM®	NP 3.2/3.3	-	-	B090-ZIA	-	Contact us
		RP 3.5/3.7/4.0/4.3	-	-	B093-ZIA	F189-ZIA	Contact us
		WP 4.5/4.6/5.0	-	-	B094-ZIA	F190-ZIA	Contact us
224	ZIACOM® External	RP	-	-	B018_01-ZIA	-	Contact us
225	ZIMMER® Eztetic™	3.1(NP)	Ti_286	-	B137	F286	S286
226	ZIMMER® Paragon™	-	-	-	-	-	-
227	ZIMMER® Swiss Plus™	4.8	Ti_109	Ti_459_NH	B060_01	F109	S109
228	ZIMMER® Spline™ A	3.25	-	-	B066_01	F110	S662
		3.75	-	-	B067_01	F111	S662
		5.0	-	-	B068_01	F112	S662
229	ZIMMER® Spline™ B	3.25	-	Ti_422_NH	B066_00	F110	S110
		3.75	-	Ti_424_NH	B067_00	F111	S110
		5.0	-	Ti_426_NH	B068_00	F112	S110
230	ZIMMER® Tapered Screw-Vent®	3.5	Ti_105	Ti_592_NH	B062_00	F157	S105
		4.5	Ti_106	Ti_594_NH	B063_00	F158	S105
		5.7	Ti_107	Ti_596_NH	B064_00	F159	S105
231	ZIMMER® Tapered Abutment SRR Multi-unit	-	-	-	M240	-	MS671
232	ZIMMER® Universal Abutment Multi-Unit	-	-	-	M200_00-ZM	MF200-ZM	Contact us

Impossible = I<sup>3</sup> m possible.



**Aligned with your vision.**



*WhatsApp*

## **Hoil Dental**

### **Address**

1 Lyon Road  
Walton-On-Thamse, KT12 3PU  
United Kingdom

### **Contact**

+44(0) 207 383 5571  
info@hoildental.com

### **Website**

www.hoildental.com  
www.hoildentalshop.com

